



Yes, you can.[®]

rea  **Azalea[®] Tall**

No limits on comfort for taller users

Rea™ **Azalea[®] Tall** is specially designed to meet the needs of tall users who require a Tilt in Space wheelchair with a longer seat support. Adapted from the Rea™ **Azalea[®]**, the Rea™ **Azalea[®] Tall** boasts all the advantages of a reliable, tilting wheelchair and offers a unique weight-shifting mechanism. The Rea™ **Azalea[®] Tall**: at the top of its class for comfort and stability.



➤ Real compatibility

The Rea™ **Azalea® Tall** has all the features of the Rea™ **Azalea®** such as Tilt in Space weight-shifting for a really stable tilt. All accessories such as rear wheels for self propelling or smaller wheels for transport and a different seat system also fit the Rea™ **Azalea® Tall**, which is in fact, simply an adaptor kit for the standard Rea™ **Azalea®** model.



➤ Easy maneuverability

Attendants will appreciate the adjustable push bar, which offers good ergonomic grip and can be adopted to fit different personal length's of the attendant.

Features



Longer seat plate

An extra 5 cm of seat depth is designed to support longer legs.



Chassis length

An extra 5 cm of height allows for longer leg-rests.



Seat height

An extra 5 cm in length enables greater forward stability.

Technical data

39-45/44-50/49-55 cm (39-49/44-54/49-59 cm)	48 - 55 cm	50 cm	50 - 80 cm (without seat cushion)	24 - 36 cm	33 - 50 cm	-1°-20° (-1°-25°)	0°-30°
SW + 22 cm (SW + 23 cm)	100 - 130 cm	95 - 107 cm	34 kg (SW 44 cm)	max. 135 kg	Transport weight 18,7 kg	Crashtested ISO 7176-19*	

Frame colour



Pearl grey

Upholstery



Grey Dartex TR23



Grey plush TR 32

* Our wheelchairs comply above ISO standards and with ISO 10542, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation. Rea, Rea design and DSS (Dual Stability System) design are registered trademarks of Invacare® International.



INVACARE® Rea AB
Box 200 · Växjövägen 303
SE-343 75 DIÖ · Sweden
Tel: +46 476 535 00 · Fax + 46 476 535 99
www.invacare-rea.com

