



Yes, you can.®

rea  **Azalea® Max**



A tailor-made solution for heavier users

Introducing a new passive (with tilt and recline) manual wheelchair for heavier users. The Rea™ **Azalea® Max** is designed to meet the specific needs of disabled users with extra weight. The **Azalea® Max** is a niche product and a complementary addition to the existing **Azalea®** family from Invacare®. Thanks to the introduction of **Azalea® Max**, dealers seeking a truly complete range of passive wheelchairs need look no further.



Continued Compatibility

The Rea™ **Azalea® Max** offers the same features as the rest of the Rea™ **Azalea®** family such as tilt-in-space with weight-shifting and DSS (Dual Stability System) mechanism. Stability in this wheelchair originates from body length compensation in the back and leg rests, beneficial for users with limited movement and that have difficulty changing position.



High ergonomic push handle

The ergonomically shaped push handles are positioned close together for easier handling. An alternate version with higher handles offers a different position for the hands and is available as an option.

Options / Features



Two-way backrest

The lower and upper parts of the shaped backrest are interchangeable, giving users the option of more space in the seat or chest area.



Footplate converter

Divided into three parts, the footplate comes with a removable mid section, to facilitate transfers.



Lower-leg support

The fixed calf plate adjusts and swivels to provide optimum support for each leg.

Technical data

Azalea® Max	55; 61 + 5cm	53 - 60	40 / 45	50 + 8 cm	32 - 42	33 - 50	1° - 20°
Azalea® Max Assist	55; 61 + 10 cm	53 - 60	40 / 45	50 + 8 cm	32 - 42	33 - 50	1° - 20°
Azalea® Max	0° - 30°	SW + 25 cm	102 - 124	105 - 112	54 kg	33 kg	max. 180 kg
Azalea® Max Assist	0° - 30°	SW + 22 cm	102 - 124	105 - 112	54 kg	37 kg	max. 180 kg

For other Rea™ **Azalea®** accessories, please consult the Rea™ **Azalea®** brochure. The Rea™ **Azalea® Max** complies with ISO 7176-19 specifications. Note that special conditions apply regarding accessories. For detailed information regarding vehicular transportation and for more pre-sale information for this wheelchair, please consult the following pages on our website: www.invacare-rea.com under Rea wheelchairs/Technical downloads.

Frame colour



Pearl grey

Upholstery



Grey Dartex TR23



Grey plush TR 32

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation. Rea, Rea design and DSS (Dual Stability System) design are registered trademarks of Invacare® International.



INVACARE® Rea AB
Box 200 · Växjövägen 303
SE-343 75 DIÖ · Sweden
Tel: +46 476 535 00 · Fax + 46 476 535 99
www.invacare-rea.com

